

Delipidation Protocol

This protocol is used to remove lipid content from fixed brain tissue via dehydration, to facilitate later section-based or whole-mount tissue staining and clearing.

Note: The optimal duration for each step depends on the sample. Example durations and buffer volumes are provided for single hemispheres and whole adult mouse brains.

Guidelines:

- For brain hemisphere: use 2mL Eppendorf tube and ~1.4-1.6mL solution for each wash.
- For whole brain: use 5mL Eppendorf tube and ~4.5mL solution for each wash.
- **Optional Tip:** During 100% MeOH and 100% DCM wash steps, move the tissue sample into an extra tube to minimize solution carry over and prevent sample from drying in the air.

IMPORTANT SAFETY NOTICE:

- 1) MeOH is volatile and flammable. Handle it with care and follow safety regulations.
- 2) DCM is volatile and an irritant. Use it inside a non-recirculating fume hood. Wear double nitrile gloves to handle the extraction process. Discard gloves immediately when contaminated with DCM.

AdipoClear+ Delipidation for Sectioning

(A) Hemisphere	(B) Whole Brain	Explanation
Total Length: 13 hrs	Total Length: 3 days	
1. Dehydrate the sample by washing up gradient of MeOH in B1n (0%, 20%, 40%, 60%, and 80% MeOH): Shake gently at RT, 30min for each step.	1. Dehydrate the sample by washing up gradient of MeOH in B1n (0%, 20%, 40%, 60%, and 80% MeOH): Shake gently at RT, 1hr for each step.	To remove relatively soluble lipids in cellular and organelle membranes.
2. Wash in 100% MeOH: shake at RT, 30min x2.	2. Wash in 100% MeOH: shake at RT, 1hr x2.	To remove hydrophobic lipids in myelin and lipid storage.
3. Wash in 100% DCM: shake at RT, 1hr.	3. Wash in 100% DCM: shake at RT, o/n.	
4. Wash in 100% MeOH: shake at RT, 30min, 45min, 1hr.	4. Wash in 100% MeOH: shake at RT, 1hr, 1.5hr, 2hr.	
5. Rehydrate with reversed MeOH/B1n gradient (from 80% to 20%): shake at RT, 30min each.	5. Rehydrate with reversed MeOH/B1n gradient (from 80% to 20%): shake at RT, 1hr each.	To rehydrate back to conditions required for immunostaining.
6. Wash in B1n buffer: shake at RT, 2hr.	7. Wash in B1n buffer: shake at RT, o/n.	

7. Wash in PBS w/ NaN_3 : shake at 37°C, 30min, 1hr, o/n.	9. Wash in PBS w/ NaN_3 : shake at 37°C, 1hr, 2hr, 4hr, o/n.	
Samples can then be stored in PBS w/ NaN_3 at 4°C for a few weeks. Fill the tube to top with media to reduce air oxidation.		