

## mAb3D Buffer Recipes

### Contents:

- I. **Suggested Reagents**
- II. **Buffer Recipes**
  - i. **Sample Collection**
    - ❖ 20mg/ml Heparin (stock)
    - ❖ 10x PBS w/ heparin (stock)
    - ❖ 4% PFA
    - ❖ B1n
  - ii. **Delipidation**
    - ❖ 10x PBS (stock)
    - ❖ 1x PBS w/ NaN<sub>3</sub>
    - ❖ MeOH/B1n gradients
    - ❖ PTxwH
  - iii. **Thick-Section Based Antibody Screening**
    - ❖ 1x PBS w/ NaN<sub>3</sub> (high conc.)
    - ❖ 20mM PB
    - ❖ AC3 mounting media

### I. Suggested Reagents

<i>Reagent</i>	<i>Reference</i>
<b>Agarose</b>	Invitrogen 16500500
<b>Dichloromethane (DCM)</b>	Sigma 270997-12X100ML
<b>Glycine</b>	Sigma G7126-500G
<b>H<sub>2</sub>O</b>	MilliQ double distilled water
<b>Heparin</b>	Sigma H3393-50KU
<b>Iohexol</b>	Nycodenz®, Axis-Shield, Cosmo Bio USA AXS-1002424 or Histodenz™, Sigma D2158-100G
<b>Methanol (MeOH)</b>	Fisher A412SK-4
<b>Paraformaldehyde 16% (PFA)</b>	EMS 15710-S
<b>Phosphate buffered saline (PBS) 10X</b>	Ambion AM9624, or prepare solution per below recipe
<b>Potassium chloride (KCl)</b>	Sigma P5405-250G
<b>Potassium phosphate monobasic (KH<sub>2</sub>PO<sub>4</sub>)</b>	Fisher P285-500
<b>Sodium azide (NaN<sub>3</sub>)</b>	Sigma 58032-100G
<b>Sodium chloride (NaCl)</b>	Fisher S640-3
<b>10N Sodium hydroxide (NaOH)</b>	Fisher SS255-1

<b>Sodium phosphate dibasic dodecahydrate (Na<sub>2</sub>HPO<sub>4</sub>·12H<sub>2</sub>O)</b>	Sigma 71650
<b>Sodium phosphate monobasic dihydrate (NaH<sub>2</sub>PO<sub>4</sub>·2H<sub>2</sub>O)</b>	Sigma 71500
<b>2,2'-Thiodiethanol (TDE)</b>	Sigma 166782-500G
<b>Triton X-100</b>	Sigma X100-500ML
<b>Tween-20</b>	Sigma P9416-100ML

## II. Buffer Recipes

### Sample Collection

#### 20mg/mL Heparin (stock): 14 mL

- Heparin salt.....280 mg
- distilled H<sub>2</sub>O.....14 mL

Add 14 mL DI H<sub>2</sub>O to 15 mL conical tube. Weigh, add heparin. Centrifuge.

#### 10x PBS w/ heparin (stock): 1 L

- Na<sub>2</sub>HPO<sub>4</sub>·12H<sub>2</sub>O.....34.2 g
- NaCl.....80.0 g
- KCl.....2.0 g
- KH<sub>2</sub>PO<sub>4</sub> .....2.0 g
- 20mg/mL Heparin.....5.0 mL
- distilled H<sub>2</sub>O

Bring to 1 L with DI H<sub>2</sub>O. Mix at RT until fully dissolved, filter; **store at RT**. To use for perfusion, dilute by factor of ten to yield 1xPBS w/ 10µg/mL heparin. Prepare buffer without heparin for tissue storage media.

#### 4% PFA: 200 mL

- 10x PBS w/ heparin.....20 mL
- PFA 16% ..... 20 mL
- distilled H<sub>2</sub>O.....130 uL

Prepare fresh each time; **keep on ice** when not in use.

#### B1n: 500 mL

- distilled H<sub>2</sub>O.....500 mL
- Glycine.....10 g
- Triton X-100.....500 uL
- NaOH, 10N.....50 uL
- 20% NaN<sub>3</sub>.....200 uL

Add water first with stir bar to prevent clumping. Add components while stirring gently. Triton is thick and sticky, need to pipette up and down; can drop pipette tip into flask to allow rest to diffuse out. Mix at RT until fully dissolved (~1hr); **store at RT**.

### Delipidation

#### **Stock solutions:**

**10x PBS: 1 L**

- Na<sub>2</sub>HPO<sub>4</sub>·12H<sub>2</sub>O.....34.2 g
- NaCl.....80.0 g
- KCl.....2.0 g
- KH<sub>2</sub>PO<sub>4</sub>.....2.0 g
- distilled H<sub>2</sub>O

Bring to 1 L with DI H<sub>2</sub>O. Mix at RT until fully dissolved; **keep at RT**.

**1x PBS w/0.02% NaN<sub>3</sub>: 1 L**

- 10x PBS.....100 mL
- 20% NaN<sub>3</sub> .....1 mL
- distilled H<sub>2</sub>O

Fill with distilled water to 1L. Mix at RT; **keep at RT**.

**Solutions for Delipidation:****MeOH gradients: 50 mL**

	20%	40%	60%	80%
MeOH	10 mL	20 mL	30 mL	40mL
B1n	40 mL	30 mL	20 mL	10 mL

Mix and cool down before use; **keep at RT**.

**Note:** MeOH/B1n series may have salt precipitation at high concentration. Use the top solution.

**PTxwH: 1L**

- distilled H<sub>2</sub>O.....900 mL
- 10x PBS.....100 mL
- Triton X-100.....1 mL
- Tween-20.....500 uL
- Heparin (20mg/ml)..... 100 uL
- 20% NaN<sub>3</sub> .....1 mL

Mix at RT until fully dissolved; **keep at RT**.

**Section-based Antibody Screening****1x PBS w/0.2% NaN<sub>3</sub> (high concentration): 50 mL**

- 10x PBS.....5.0 mL
- 20% NaN<sub>3</sub> .....0.5 mL
- distilled H<sub>2</sub>O.....45.0 mL

Mix at RT; **keep at RT**. **Use:** to initially fill well plates for sections.

**20 mM PB: 1L**

- Na<sub>2</sub>HPO<sub>4</sub>·12H<sub>2</sub>O...(sodium phosphate dibasic dodecahydrate).....5.73 g
- NaH<sub>2</sub>PO<sub>4</sub>·2H<sub>2</sub>O.....(sodium phosphate monobasic dihydrate).....0.624 g
- distilled H<sub>2</sub>O

Bring to 1 L with DI H<sub>2</sub>O. Mix at RT until fully dissolved; **keep at RT**.

**ACB high RI mounting media:** in 50ml conical tube

- 20 mM PB.....8.0 mL
- 2,2'-Thiodiethanol (TDE).....22.4 mL
- Iohexol (Histodenz or Nycodenz).....24.0 g

Mix the PB buffer and TDE first to prevent clumping. Shake at RT, 200 rpm o/n, until dissolved. Spin down at 5,000 rpm for 10 min to remove any precipitate. Aliquot for use at **RT**.